**BioPoly<sup>®</sup>**

# Partial Resurfacing Patella Implant

**Knieerneuerungs-Teilimplantat**  
**Implant de resurfaçage partiel de la rotule**  
**Implante rotuliano para reconstrucción parcial**  
**Impianto di ricostruzione parziale superficiale della rotula**  
**Patellaimplantat til delvis overfladeudskiftning**  
**Knieschijfimplantaat voor gedeeltelijke resurfacing**  
**Polvilumpion osittaisessa uudelleenpinnoittamisessa käytettävä implantti**  
**Implant do częściowej protezoplastyki powierzchniowej rzepki**  
**Implante de Revestimento Parcial de Rótula**  
**Partiellt knäledsimplantat**  
**Kısmi Yüzey Değişirme Patella İmplantı**  
**Implantát pro částečnou povrchovou náhradu česky**

Instructions for Use  
Gebrauchsanweisung  
Mode d'emploi  
Instrucciones de uso  
Istruzioni per l'uso  
Brugsvejledning  
Gebruiksaanwijzing  
Käyttöohjeet  
Instrukcja użycia  
Instruções de Utilização  
Bruksanvisning  
Kullanım Talimatları  
Návod k použití

Manufacturer

BioPoly LLC  
7136 Gettysburg Pike  
Fort Wayne, IN 46804  
USA  
+1-260-999-6135

STERILE EO

Sterilized Using  
Ethylene Oxide

Single Use

Attention, See  
Instructions For Use

125-3004-01 rev 1.5  
Date: 2016-12

## ENGLISH

### DEVICE DESCRIPTION

The BioPoly® RS Partial Resurfacing Patella Implant is a long term, surgically invasive device for replacement of focal osteochondral defects in the patella less than 3.1cm<sup>2</sup>.

Caution: For Cemented Use Only

### MATERIALS

The implant is made of BioPoly® material (ultra-high molecular weight polyethylene and cross-linked hyaluronan).

### INTENDED USE

The BioPoly® RS Partial Resurfacing Patella Implant is intended for the replacement of symptomatic abnormal or severely abnormal (ICRS Grade 2, 3, or 4) chondral or osteochondral focal lesions located in the medial and/or lateral facet of the patella in patients over 21 years of age. Lesion size may not exceed 3.1cm<sup>2</sup> and must be circumscribed by a 1.5 or 2.0cm circle of normal or nearly normal (ICRS Grade 0 or 1) cartilage, with an overall cartilage depth less than 4mm from the articulating surface. Subchondral bone quality must be deemed sufficient to secure the implant with bone cement.

### CONTRAINDICATIONS

1. Body mass index  $\geq 35$
2. Generalized degenerative or autoimmune arthritis
3. Gout
4. Uncorrected malalignment of the patella
5. Uncorrected ligamentous instability
6. Kissing lesion on femur
7. More than one implant required to accommodate lesion
8. Allergy to ultra-high molecular weight polyethylene (UHMWPE), or hyaluronic acid (HA)
9. Use with opposing articulating femoral components

### WARNINGS

1. Prior to use, the surgeon should thoroughly read and understand all aspects of the surgical technique.
2. This device is provided STERILE as a single use product. If sterile barrier is broken or the packaging otherwise damaged, the device may not be administered. Re-sterilization of the device is not permitted.
3. During removal of the device from the packaging, special care should be taken to preserve the articulating surface finish.
4. In case of damage to the implant, including scratches or indentations on the articulating surface, bending of the fixation stem, or any other disfigurement, the device may not be administered.
5. Any alteration or modification to the device prior to the surgical implantation is prohibited.
6. Excessive insertion force can damage the articulating surface of the implant.
7. Do not reuse device. Risks of reuse include: damage to implant, loss of performance, and infection.
8. Implant articulation with non cartilage or abnormal anatomic surface can damage the implant.

### ADVERSE EFFECTS

Infections, allergies, or other reactions to implant materials.

### PACKAGING AND STERILIZATION

The BioPoly® RS Partial Resurfacing Patella implant is provided sterile and for single use.

1. Do not use if the sterile packaging has been breached or damaged
2. Do not attempt to re-sterilize the implant
3. Do not use if the expiration date has elapsed

### DIRECTIONS FOR USE

The BioPoly® Instrumentation must be used to ensure proper fit and alignment of the implant. Refer to the appropriate surgical technique manual for specific instructions on the use of the instruments and implantation of the prosthesis

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