BioPoly[®]

Partial Resurfacing Patella / Central Ridge Patella Implant

Teilerneuerung der Patella / Zentraler Kamm des Patella-Implantats
Implant de Resurfaçage Partiel de la Rotule / Arête Centrale de la Rotule
Implante Rotuliano Patelar / Cresta Central para Reconstrucción Parcial
Impianto di Ricostruzione Parziale Superficiale della Rotula / Cresta Centrale
Patellaimplantat / Central Ryg Patellaimplantat til delvis overfladeudskiftning
Knieschijf / Centrale Nokimplantaat voor Gedeeltelijke Resurfacing
Implantti Patella / Keskihajanteen Osittaiseksi Uudelleen Pinnoittamiseksi
Częściowej Implant nawierzchniowy Rzepki / Grzbietu Centralnego
Implante de Revestimento Parcial de Rótula / Cresta Central de Rótula
Partiellt Partiell återupptäckande Patella / Central Ryggpatellaimplantat
Kısmi Yüzey Değiştirme Dizkapağı / Patellanın Orta Sırtı İmplantı
Implantát pro částečnou Výměnu Povrchu Hřebenu Patelly / Patelly
Delvis Oppblomstring Patella / Sentral ås på Patellaimplantat

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Manufacturer
BioPoly LLC
7136 Gettysburg Pike
Fort Wayne, IN 46804
USA
+1-260-999-6135

European Authorized
Representative

HealthLink Europe
De Tweeling 20-22
5215 MC 's-
Hertogenbosch
The Netherlands

Sterilized Using
Ethylene Oxide

Single Use

Attention, See
Instructions For Use

125-3004-01 rev 1.7
Date: 2019-08

DEVICE DESCRIPTION

The BioPoly® RS Partial Resurfacing Patella / Central Ridge Patella Implant is a long term, surgically invasive device for replacement of focal osteochondral defects in the patella less than 3.1cm².

Caution: For Cemented Use Only

MATERIALS

The implant is made of BioPoly® material (ultra-high molecular weight polyethylene and cross-linked hyaluronan).

INTENDED USE

The BioPoly® RS Partial Resurfacing Patella / Central Ridge Patella Implant is intended for the replacement of symptomatic abnormal or severely abnormal (ICRS Grade 2, 3, or 4) chondral or osteochondral focal lesions located in the medial facet, lateral facet, or central ridge of the patella in patients over 21 years of age. Lesion size may not exceed 3.1cm² and must be circumscribed by a 1.5 or 2.0cm circle of normal or nearly normal (ICRS Grade 0 or 1) cartilage, with an overall cartilage depth less than 4mm from the articulating surface. Subchondral bone quality must be deemed sufficient to secure the implant with bone cement.

CONTRAINDICATIONS

1. Body mass index ≥ 35
2. Generalized degenerative or autoimmune arthritis
3. Gout
4. Uncorrected malalignment of the patella
5. Uncorrected ligamentous instability
6. Kissing lesion on femur
7. More than one implant required to accommodate lesion
8. Allergy to ultra-high molecular weight polyethylene (UHMWPE), or hyaluronic acid (HA)
9. Use with opposing articulating femoral components

WARNINGS

1. Prior to use, the surgeon should thoroughly read and understand all aspects of the surgical technique.
2. This device is provided STERILE as a single use product. If sterile barrier is broken or the packaging otherwise damaged, the device may not be administered. Re-sterilization of the device is not permitted.
3. During removal of the device from the packaging, special care should be taken to preserve the articulating surface finish.
4. In case of damage to the implant, including scratches or indentations on the articulating surface, bending of the fixation stem, or any other disfigurement, the device may not be administered.
5. Any alteration or modification to the device prior to the surgical implantation is prohibited.
6. Excessive insertion force can damage the articulating surface of the implant.
7. Do not reuse device. Risks of reuse include: damage to implant, loss of performance, and infection.
8. Implant articulation with non cartilage or abnormal anatomic surface can damage the implant.
9. For Anatomic Implant, ensure proper orientation is followed per surgical technique.

ADVERSE EFFECTS

Infections, allergies, or other reactions to implant materials.

PACKAGING AND STERILIZATION

The BioPoly® RS Partial Resurfacing Patella / Central Ridge Patella implant is provided sterile and for single use.

1. Do not use if the sterile packaging has been breached or damaged
2. Do not attempt to re-sterilize the implant
3. Do not use if the expiration date has elapsed

DIRECTIONS FOR USE

Refer to surgical technique for selection of appropriate implant size and shape. The BioPoly® Instrumentation must be used to ensure proper fit and alignment of the implant. Refer to the appropriate surgical technique for specific instructions on the use of the instruments and implantation of the prosthesis

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