BioPoly[®]

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USA

+1-260-999-6135

European Authorized
Representative

HealthLink Europe

De Tweeling 20-22

5215 MC's-Hertogenbosch

The Netherlands

Sterilized
Using Ethylene
Oxide

Single
Use

Attention, See
Instructions For
Use

120-3004-01 rev 1.2
Date: 2019-07

DEVICE DESCRIPTION

The BioPoly® -RS Partial Resurfacing Shoulder Implant is a long term, surgically invasive device for replacement of focal osteochondral defects in the shoulder less than 7cm².

MATERIALS

The implant is made of BioPoly® material (ultra-high molecular weight polyethylene and crosslinked hyaluronan) and titanium alloy (Ti-6Al-4V).

INTENDED USE

The BioPoly® -RS Partial Resurfacing Shoulder Implant is intended for the replacement of symptomatic abnormal or severely abnormal (ICRS Grade 2, 3, or 4) chondral or osteochondral focal lesions located in the humeral head in patients over 21 years of age. Lesion size may not exceed 7cm² and must be circumscribed by a 25 or 30mm circle of normal or nearly normal (ICRS Grade 0 or 1) cartilage, with an overall cartilage depth less than 4mm from the articulating surface. For use without bone cement. Subchondral bone quality must be deemed sufficient to secure the implant.

CONTRAINDICATIONS

1. Generalized degenerative or autoimmune arthritis
2. Gout
3. Uncorrected chronic malalignment of the shoulder
4. Uncorrected ligamentous instability
5. Uncorrected mechanically symptomatic labral tear
6. Kissing lesion on glenoid
7. More than one implant required to accommodate lesion
8. Allergy to titanium alloy (Ti-6Al-4V), ultra-high molecular weight polyethylene (UHMWPE), or hyaluronic acid (HA)
9. Use with opposing articulating glenoid components

WARNINGS

1. Prior to use, the surgeon should thoroughly read and understand all aspects of the surgical technique.
2. This device is provided STERILE as a single use product. If sterile barrier is broken or the packaging otherwise damaged, the device may not be administered. Re-sterilization of the device is not permitted.
3. During removal of the device from the packaging, special care should be taken to preserve the articulating surface finish.
4. In case of damage to the implant, including scratches or indentations on the articulating surface, bending of the fixation stem, or any other disfigurement, the device may not be administered.
5. Any alteration or modification to the device prior to the surgical implantation is prohibited.
6. Excessive insertion force can damage the articulating surface of the implant.
7. Implant articulation with non cartilage or abnormal anatomic surface can damage the implant.
8. Do not reuse device. Risks of reuse include: damage to implant, loss of performance, and infection.
9. MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device.

ADVERSE EFFECTS

Infections, allergies, or other reactions to implant materials.

PACKAGING AND STERILIZATION

The BioPoly® Partial Resurfacing Shoulder implant is provided sterile and for single use.

1. Do not use if the sterile packaging has been breached or damaged
2. Do not attempt to re-sterilize the implant
3. Do not use if the expiration date has elapsed

DIRECTIONS FOR USE

The BioPoly® instrumentation must be used to ensure proper fit and alignment of the implant. Refer to the appropriate surgical technique manual for specific instructions on the use of the instruments and implantation of the prosthesis.

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