BioPoly®

Partial Resurfacing Knee Implant

Knieerneuerungs-Teilimplantat
Implant de resurfaçage partiel du genou
Implante de rodilla de restauración parcial
Impianto di rivestimento parziale del ginocchio
Knæimplantat til delvis rekonstruktion
Knie-implantaat voor gedeeltelijke resurfacing
Polven osittainen uudelleenpinnoitusimplantti
Implant do częściowej protezoplastyki
powierzchniowej kolana
Implante de Revestimento Parcial do Joelho
Partiella knäledsimplantat
Kısmi Yüzey Değiştirme Diz İmplantı
Implantát pro částečnou povrchovou náhradu
konce kolenního kloubu
Delvis Resurfacing Knæ Implantat
Parțial Implantul de Genunchi
Delna površina kolenskega vsadka
Čiastočný resurfacing kolenného implantátu
Делимично враћање колена имплантата

Partial Resurfacing Trochlea Implant

Trochlea erneuerungs-Teilimplantat
Implant de resurfaçage partiel de la trochlée
Implante de troclear restaurativo parcial
Impianto di rivestimento parziale del trochlea
Trochlea implantat til delvis rekonstruktion
Trochleair-implantaat voor gedeeltelijke resurfacing
Trochlea osittainen uudelleenpinnoitusimplantti
Implant do częściowej protezoplastyki Par
Implante de Revestimento Parcial do Trochlea
Partiella trochleaimplantat
Kısmi Yüzey Değiştirme trokleası İmplantları
Implantát pro částečnou výměnu povrchu trochlea
Delvis Resurfacing Trochlea Implantat
Implantare trochlea cu suprafață parțială
Delna površinska implantacija trochlea
Čiastočný resurfacing Trochlea Implant
Партиал Ресурфацинг Троцхлеа Имплант

Instructions for Use
Gebrauchsanweisung
Mode d'emploi
Instrucciones de uso
Istruzioni per l'uso
Brugsvejledning
Gebruiksaanwijzing
Käyttöohjeet
Instrukcja użycia

Instruções de Utilização
Bruksanvisning
Kullanım Talimatları
Návod k použití
Instruksjoner til bruk
Instructiuni de folosire
Navodila za uporabo
Inštrukcie na používanie
упутство за употребу

Manufacturer

BioPoly LLC
7136 Gettysburg Pike
Fort Wayne, IN 46804
USA
+1-260-999-6135

STERILE EO

Sterilized Using
Ethylene Oxide

Single Use

Attention, See
Instructions For Use

115-3004-01 rev 2.3
Date: 2018-08

ENGLISH

DEVICE DESCRIPTION

The BioPoly® RS Partial Resurfacing Knee/Trochlea Implant is a long term, surgically invasive device for replacement of focal osteochondral defects in the knee ranging from 0.8 – 3.1cm².

MATERIALS

The implant is made of BioPoly® material (ultra-high molecular weight polyethylene and crosslinked hyaluronan) and titanium alloy (Ti-6Al-4V).

INTENDED USE

The BioPoly® RS Partial Resurfacing Knee/Trochlea Implant is intended for the replacement of symptomatic abnormal or severely abnormal (ICRS Grade 2, 3, or 4) chondral or osteochondral focal lesions located in the weight bearing regions of the femoral condyle, the femoral trochlear facets, or the femoral trochlear groove in patients over 21 years of age. Lesion size may not exceed 3.1cm² and must be circumscribed by a 1.0, 1.5, or 2.0cm circle or 1.5 by 2.4cm oval (femoral condyles and trochlear facets only) of normal or nearly normal (ICRS Grade 0 or 1) cartilage, with an overall cartilage depth less than 4mm from the articulating surface. For use with or without bone cement. Subchondral bone quality must be deemed sufficient to secure the implant.

CONTRAINDICATIONS

1. Body mass index ≥ 35
2. Generalized degenerative or autoimmune arthritis
3. Gout
4. Uncorrected chronic malalignment of the knee
5. Uncorrected ligamentous instability
6. Uncorrected mechanically symptomatic meniscal tear or total meniscectomy
7. Kissing lesion on tibia or patella
8. More than one implant required to accommodate lesion
9. Allergy to titanium alloy (Ti-6Al-4V), ultra-high molecular weight polyethylene (UHMWPE), or hyaluronic acid (HA)
10. Use with opposing articulating tibial or patellar components

WARNINGS

1. Prior to use, the surgeon should thoroughly read and understand all aspects of the surgical technique.
2. This device is provided STERILE as a single use product. If sterile barrier is broken or the packaging otherwise damaged, the device may not be administered. Re-sterilization of the device is not permitted.
3. During removal of the device from the packaging, special care should be taken to preserve the articulating surface finish.
4. In case of damage to the implant, including scratches or indentations on the articulating surface, bending of the fixation stem, or any other disfigurement, the device may not be administered.
5. Any alteration or modification to the device prior to the surgical implantation is prohibited.
6. Excessive insertion force can damage the articulating surface of the implant.
7. Implant articulation with non cartilage or abnormal anatomic surface can damage the implant.
8. Do not reuse device. Risks of reuse include: damage to implant, loss of performance, and infection.
9. MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device.
10. Care should be taken to ensure that the BioPoly® RS Partial Resurfacing Knee implant is not placed posteriorly beyond the weight bearing region when using in the femoral condyle due to potential mismatch with anatomy.

ADVERSE EFFECTS

Infections, allergies, or other reactions to implant materials.

PACKAGING AND STERILIZATION

The BioPoly® Partial Resurfacing Knee implant is provided sterile and for single use.

1. Do not use if the sterile packaging has been breached or damaged
2. Do not attempt to re-sterilize the implant
3. Do not use if the expiration date has elapsed

DIRECTIONS FOR USE

The BioPoly® instrumentation must be used to ensure proper fit and alignment of the implant. Refer to the appropriate surgical technique manual for specific instructions on the use of the instruments and implantation of the prosthesis.

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