



BioPoly[®] RS

Patella System


Patella partial resurfacing device
made from a microcomposite of
Hyaluronic Acid and UHMWPE



Bio
Hyaluronic Acid



Poly
UHMWPE



BioPoly[®]


Hydrophilic polymer provides
a smooth and lubricious
surface to interact favorably
with native joint tissues and
carry anatomical loads

BioPoly RS Patella Implant




Patient presents with a painful cartilage lesion of the media and/or lateral facet of the patella.

- Arthroscopic assisted implantation is quick and straightforward.
- Joint disruption is minimized.
- Anatomy is preserved.
- Unique implant design and construction restores articulating surface.
- Patient rapidly recovers and returns to activity.



Reference No.	Description
125-5004	BioPoly 15mm Partial Resurfacing Patella Implant (Size 0)
125-5005	BioPoly 15mm Partial Resurfacing Patella Implant (Size 1)
125-5006	BioPoly 20mm Partial Resurfacing Patella Implant (Size 0)
125-5007	BioPoly 20mm Partial Resurfacing Patella Implant (Size 1)
115-5500	Instrument Kit



15 mm
Size 0

15 mm
Size 1

20 mm
Size 0

20 mm
Size 1

BioPoly® RS device is implanted utilizing a reusable instrument kit. Please contact a BioPoly, LLC representative to learn more about the benefits of the BioPoly RS Patella System. Sales@BioPolyortho.com

